

Anti-TTC7B Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A12170-1-Fluoro594

About TTC7B

Component of a complex required to localize phosphatidylinositol 4-kinase(PI4K) to the plasma membrane. The complex acts as a regulator of phosphatidylinositol 4-phosphate(PtdIns(4)P) synthesis. In the complex, plays a central role in bridging PI4KA to EFR3B and FAM126A, via direct interactions

Overview

Product Name	Anti-TTC7B Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q86TV6

Technical Details

Immunogen	E.coli-derived human TTC7B recombinant protein (Position: E275-H730).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TTC7B Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.