

Anti-EEFSEC Antibody Picoband® Biotin Conjugated

Catalog Number: A12172-2-Biotin

About EEFSEC

Enables GTPase activity. Involved in selenocysteine incorporation. Predicted to be located in cytoplasm and nucleus. Predicted to be part of ribonucleoprotein complex.

Overview

Product Name	Anti-EEFSEC Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P57772

Technical Details

Immunogen	E.coli-derived human EEFSEC recombinant protein (Position: P216-P596).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-EEFSEC Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.