

## Anti-ABT1 Antibody (C-term)

Catalog Number: A12254

### About ABT1

Could be a novel TATA-binding protein (TBP) which can function as a basal transcription activator. Can act as a regulator of basal transcription for class II genes (By similarity).


### Overview

Product Name	Anti-ABT1 Antibody (C-term)
Reactive Species	Human
Description	Boster Bio Anti-ABT1 Antibody (C-term) (Catalog # A12254). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9ULW3

### Technical Details

Immunogen	This ABT1 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 259-292 amino acids from the C-terminal region of human ABT1.
Predicted Reactive Species	Bovine, Mouse, Rat
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Suggested Dilutions	WB: 1:1000

## Anti-ABT1 Antibody (C-term) (A12254) Images



Western blot analysis of lysate from HepG2 cell line, using ABT1 Antibody (C-term). A12254 was diluted at 1:1000. A goat anti-rabbit IgG H&L (HRP) at 1:10000 dilution was used as the secondary antibody. Lysate at 20ug.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ABT1 Antibody (C-term)

For Research Use Only. Not for use in diagnostic procedures.