

## Anti-PARP15 Antibody Picoband® FITC Conjugated

Catalog Number: A12266-1-FITC

### About PARP15

Enables NAD<sup>+</sup> binding activity; pentosyltransferase activity; and transcription corepressor activity. Involved in negative regulation of transcription by RNA polymerase II and protein poly-ADP-ribosylation. Predicted to be active in cytoplasm and nucleus. Biomarker of atherosclerosis.

### Overview

Product Name	Anti-PARP15 Antibody Picoband® FITC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q460N3

### Technical Details

Immunogen	E.coli-derived human PARP15 recombinant protein (Position: N47-R555). Human PARP15 shares 46.6% amino acid (aa) sequence identity with mouse PARP15.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-PARP15 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.