

Anti-PARP15 Antibody Picoband® HRP Conjugated

Catalog Number: A12266-1-HRP

About PARP15

Enables NAD⁺ binding activity; pentosyltransferase activity; and transcription corepressor activity. Involved in negative regulation of transcription by RNA polymerase II and protein poly-ADP-ribosylation. Predicted to be active in cytoplasm and nucleus. Biomarker of atherosclerosis.

Overview

Product Name	Anti-PARP15 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q460N3

Technical Details

Immunogen	E.coli-derived human PARP15 recombinant protein (Position: N47-R555). Human PARP15 shares 46.6% amino acid (aa) sequence identity with mouse PARP15.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PARP15 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.