

Anti-DDI2 Antibody Picoband® FITC Conjugated

Catalog Number: A12277-1-FITC

About DDI2

Enables aspartic-type endopeptidase activity; identical protein binding activity; and ubiquitin binding activity. Involved in several processes, including cellular response to hydroxyurea; proteolysis; and regulation of DNA stability. Located in cytosol and nucleoplasm.

Overview

Product Name	Anti-DDI2 Antibody Picoband® FITC Conjugated
Reactive Species	Human, Mouse
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IP, IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q5TDH0

Technical Details

Immunogen	E.coli-derived human DDI2 recombinant protein (Position: Q46-P399).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-DDI2 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.