

Anti-HRP-3 HDGFL3 Antibody

Catalog Number: A12283-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-HRP-3 HDGFL3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-HRP-3 HDGFL3 Antibody catalog # A12283-1. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9Y3E1

Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human HDGFRP3. AA range:131-180
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000

	IHC: 1:100-1:300 ELISA 1:20000 IF 1:50-200 Specific for endogenous levels of the ~62 kDa COBRA1 protein.
--	---

Anti-HRP-3 HDGFL3 Antibody (A12283-1) Images



Western Blot analysis of NIH-3T3, PC12 cells using HRP-3 Polyclonal Antibody. Secondary antibody was diluted at 1:20000

Anti-HRP-3 HDGFL3 Antibody

For Research Use Only. Not for use in diagnostic procedures.