

## Anti-NCLN Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A12284-2-Fluoro647

### About NCLN

Nicalin is a protein that in humans is encoded by the NCLN gene. Enables ribosome binding activity. Involved in protein stabilization; regulation of protein complex stability; and regulation of protein-containing complex assembly. Located in endoplasmic reticulum membrane. Part of protein-containing complex.

### Overview

Product Name	Anti-NCLN Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q969V3

### Technical Details

Immunogen	E.coli-derived human NCLN recombinant protein (Position: Q59-Q563). Human NCLN shares 92.9% and 92.5% amino acid (aa) sequence identity with mouse and rat NCLN, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-NCLN Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.