

Anti-HDDC3 Antibody Picoband® Cy3 Conjugated

Catalog Number: A12287-1-Cy3

About HDDC3

HDDC3, also known as Histone Deacetylase Domain Containing 3, is a gene that encodes a protein with a potential histone deacetylase activity. Proteins in this family are involved in the modification of histones, which can affect gene expression and cellular processes. Although the precise functions of HDDC3 are not fully understood, it is likely to play a role in regulating gene expression and cellular differentiation.

Overview

Product Name	Anti-HDDC3 Antibody Picoband® Cy3 Conjugated
Reactive Species	Human, Mouse
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8N4P3

Technical Details

Immunogen	E.coli-derived human HDDC3 recombinant protein (Position: K25-I179). Human HDDC3 shares 95.5% amino acid (aa) sequence identity with mouse HDDC3.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HDDC3 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.