

Anti-NDF6 NEUROD6 Antibody

Catalog Number: A12294


Overview

Product Name	Anti-NDF6 NEUROD6 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-NDF6 NEUROD6 Antibody catalog # A12294. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q96NK8


Technical Details

Immunogen	Synthesized peptide derived from human protein at AA range: 230-310
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	NDF6 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 ELISA 1:5000-20000

Anti-NDF6 NEUROD6 Antibody (A12294) Images



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4°over night



Western Blot analysis of Hela lysis, using primary antibody at 1:1000 dilution. Secondary antibody was diluted at 1:10000

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NDF6 NEUROD6 Antibody

For Research Use Only. Not for use in diagnostic procedures.