

Anti-UCHL5IP/HAUS7 Antibody Picoband® Biotin Conjugated

Catalog Number: A12338-1-Biotin

About HAUS7

This gene encodes a subunit of the augmin complex, which regulates centrosome and mitotic spindle integrity, and is necessary for the completion of cytokinesis. The encoded protein was identified by interaction with ubiquitin C-terminal hydrolase 37. Alternative splicing results in multiple transcript variants.

Overview

Product Name	Anti-UCHL5IP/HAUS7 Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q99871

Technical Details

Immunogen	E.coli-derived human UCHL5IP/HAUS7 recombinant protein (Position: D18-Q333). Human HAUS7 shares 69% amino acid (aa) sequence identity with mouse HAUS7.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-UCHL5IP/HAUS7 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.