

## Anti-Ras-related protein Rab-9B RAB9B Antibody

Catalog Number: A12384


### Overview

|                      |   |
|----------------------|---|
| Product Name         | Anti-Ras-related protein Rab-9B RAB9B Antibody  |
| Reactive Species     | Human, Mouse  |
| Description          | Boster Bio Anti-Ras-related protein Rab-9B RAB9B Antibody catalog # A12384. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse. |
| Application          | ELISA, WB   |
| Clonality            | Polyclonal  |
| Formulation          | Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.   |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.            |
| Host                 | Rabbit  |
| Uniprot ID           | Q9NP90  |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | Synthesized peptide derived from human protein at AA range: 80-160  |
| Isotype             | IgG   |
| Form                | Liquid  |
| Concentration       | 1 mg/ml   |
| Purification        | RAB9B antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Suggested Dilutions | WB 1:500-2000<br>ELISA 1:5000-20000   |

## Anti-Ras-related protein Rab-9B RAB9B Antibody (A12384) Images



Western blot analysis of lysates from Jurkat cells, primary antibody was diluted at 1:1000, 4°over night

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Ras-related protein Rab-9B RAB9B Antibody

For Research Use Only. Not for use in diagnostic procedures.