

Anti-Cadherin-9 (P222) CDH9 Antibody

Catalog Number: A12410-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-Cadherin-9 (P222) CDH9 Antibody
Reactive Species	Rat
Description	Boster Bio Anti-Cadherin-9 (P222) CDH9 Antibody catalog # A12410-1. Tested in IHC applications. This antibody reacts with Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9ULB4

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 190-241 of Human Cadherin-9.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

Anti-Cadherin-9 (P222) CDH9 Antibody (A12410-1) Images



Immunohistochemistry (IHC) analyzes of Cadherin-9 (P222) pAb in paraffin-embedded human brain tissue.

[Anti-Cadherin-9 \(P222\) CDH9 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.