

Anti-ARHGAP25 Antibody

Catalog Number: A12463-1

About ARHGAP25

ARHGAP25 (Rho GTPase-activating protein 25) also known as KAIA0053, belongs to a family of Rho GTPase-modulating proteins which are implicated in actin remodeling, cell polarity, and cell migration. ARHGAP25 has a key role in the pathomechanism of autoantibody-induced arthritis in which it regulates inflammation.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-ARHGAP25 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ARHGAP25 Antibody catalog # A12463-1. Tested in WB, IHC, ICC/IF, IP, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P42331

Technical Details

Immunogen	E.coli-derived human ARHGAP25 recombinant protein (Position: L3-E344).
Form	Liquid

Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 ImmunoPrecipitation, 1:250-300 ELISA, 1:100-1000

Anti-ARHGAP25 Antibody

For Research Use Only. Not for use in diagnostic procedures.