

Anti-Histone H3 (phospho-S10) HIST1H3A Antibody

Catalog Number: A12477S10


Overview

Product Name	Anti-Histone H3 (phospho-S10) HIST1H3A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Histone H3 (phospho-S10) HIST1H3A Antibody catalog # A12477S10. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P68431/Q71DI3/P84243


Technical Details

Immunogen	Synthetic phosphopeptide derived from human Histone H3 around the phosphorylation site of Serine 10.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IF: 1:50-1:200

Anti-Histone H3 (phospho-S10) HIST1H3A Antibody (A12477S10) Images



Western blot (WB) analysis of p-Histone H3 (S10) pAb at 1:1000 dilution Lane1:HEK293T whole cell lysate(40ug) Lane2:MCF-7 whole cell lysate(10ug) Lane3:PC3 whole cell lysate(40ug) Lane4:PC12 whole cell lysate(40ug) Lane5:AML-12 whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of p-Histone H3 (S10) pAb in paraffin-embedded human skeletal muscle tissue.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Histone H3 (phospho-S10) HIST1H3A Antibody

For Research Use Only. Not for use in diagnostic procedures.