

Anti-Histone H3 (Phospho-S10) HIST1H3A Antibody

Catalog Number: A12477S10-1


Overview

Product Name	Anti-Histone H3 (Phospho-S10) HIST1H3A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Histone H3 (Phospho-S10) HIST1H3A Antibody catalog # A12477S10-1. Tested in WB,ICC/IF,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P68431


Technical Details

Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Ser10 of human Histone H3.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:5,000 ICC: 1:50-1:200 IHC: 1:50-1:200

Anti-Histone H3 (Phospho-S10) HIST1H3A Antibody (A12477S10-1) Images



Western blot analysis of Phospho-Histone H3(S10) on HeLa cell lysates using anti-Phospho-Histone H3(S10) antibody at 1/1,000 dilution. Positive control: Lane 1: Untreated HeLa cell lysate Lane 2: HeLa cell lysate treated with calyculin A Lane 3: HeLa cell lysate treated with calyculin A and alkaline phosphatase



ICC staining Phospho-Histone H3(S10) in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Histone H3 (Phospho-S10) HIST1H3A Antibody

For Research Use Only. Not for use in diagnostic procedures.