

## Anti-PAF65 Alpha (Q53) TAF6L Antibody

Catalog Number: A12489-1


### Overview

Product Name	Anti-PAF65 Alpha (Q53) TAF6L Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-PAF65 Alpha (Q53) TAF6L Antibody catalog # A12489-1. Tested in WB,IHC applications. This antibody reacts with Human,Mouse.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9Y6J9


### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 20-70 of Human PAF65alpha.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

## Anti-PAF65 Alpha (Q53) TAF6L Antibody (A12489-1) Images



Western blot (WB) analysis of PAF65alpha (Q53) Antibody at 1:1000 dilution Lane1:HepG2 whole cell lysate(10ug) Lane2:HEK293T whole cell lysate(10ug) Lane3:A549 whole cell lysate(10ug) Lane4:EC9706 whole cell lysate(10ug) Lane5:The Lung tissue lysate of Mouse(40ug)



Immunohistochemistry (IHC) analyzes of PAF65alpha (Q53) pAb in paraffin-embedded human tonsil cancer tissue at 1:50.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PAF65 Alpha (Q53) TAF6L Antibody

For Research Use Only. Not for use in diagnostic procedures.