

# Anti-PAF65 Alpha (Q53) TAF6L Antibody

Catalog Number: A12489-1


## Overview

|                      |                                                                                                                                            |
|----------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-PAF65 Alpha (Q53) TAF6L Antibody                                                                                                      |
| Reactive Species     | Human, Mouse                                                                                                                               |
| Description          | Boster Bio Anti-PAF65 Alpha (Q53) TAF6L Antibody catalog # A12489-1. Tested in WB,IHC applications. This antibody reacts with Human,Mouse. |
| Application          | IHC, WB                                                                                                                                    |
| Clonality            | Polyclonal                                                                                                                                 |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2                                                                     |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host                 | Rabbit                                                                                                                                     |
| Uniprot ID           | Q9Y6J9                                                                                                                                     |


## Technical Details

|                     |                                                                                                                                                             |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen           | Synthetic peptide, corresponding to amino acids 20-70 of Human PAF65alpha.                                                                                  |
| Isotype             | IgG                                                                                                                                                         |
| Form                | Liquid                                                                                                                                                      |
| Concentration       | 1 mg/ml                                                                                                                                                     |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:1000<br>IHC: 1:50-1:200                                                                                                                         |

## Anti-PAF65 Alpha (Q53) TAF6L Antibody (A12489-1) Images



Western blot (WB) analysis of PAF65alpha (Q53) Antibody at 1:1000 dilution Lane1:HepG2 whole cell lysate(10ug) Lane2:HEK293T whole cell lysate(10ug) Lane3:A549 whole cell lysate(10ug) Lane4:EC9706 whole cell lysate(10ug) Lane5:The Lung tissue lysate of Mouse(40ug)



Immunohistochemistry (IHC) analyzes of PAF65alpha (Q53) pAb in paraffin-embedded human tonsil cancer tissue at 1:50.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PAF65 Alpha (Q53) TAF6L Antibody

For Research Use Only. Not for use in diagnostic procedures.