

## Anti-Twistnb Antibody

Catalog Number: A12525

### About TWISTNB

This antibody is suitable for Cancer, Immunology and Nuclear Signaling research. TMBIM1 (Transmembrane BAX inhibitor motif-containing protein 1) is a member of the 7 TMS (7 transmembrane domains) family of receptors known to mediate the activation of various transcription factors. TMBIM1 is identified as a novel modulator of NF- $\kappa$ B activation. A drosophila homolog of this protein (dNMDA1 with 40% homology to hTMBIM1) is reported to be unregulated during aging and oxidative stress.

### Overview

Product Name	Anti-Twistnb Antibody
Reactive Species	Human
Description	Boster Bio Anti-Twistnb Antibody (Catalog# A12525). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	100 $\mu$ g in 200 $\mu$ l PBS containing 0.05% stabilizing protein and 0.05% sodium azide This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store the antibody at 4°C, stable for 6 months. For long-term storage, store at -20°C. Avoid repeated freeze and thaw cycles.
Host	Rabbit
Uniprot ID	Q3B726


### Technical Details

Immunogen	A partial length recombinant Twistnb protein (amino acids 28-257)
Predicted Reactive Species	Chimpanzee
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Protein A Chromatography

Suggested Dilutions

Western blot: 1-2  $\mu\text{g/ml}$

## Anti-Twistnb Antibody (A12525) Images



Expression analysis of Twistnb. Anti-Twistnb antibody was used at 1 &mu

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-Twistnb Antibody

For Research Use Only. Not for use in diagnostic procedures.