

Anti-LPPR4 PLPPR4 Antibody

Catalog Number: A12545

Overview

Product Name	Anti-LPPR4 PLPPR4 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-LPPR4 PLPPR4 Antibody catalog # A12545. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q7Z2D5

Technical Details

Immunogen	Synthesized peptide derived from part region of human protein.AA55-80
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	LPPR4 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 ELISA 1:5000-20000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LPPR4 PLPPR4 Antibody

For Research Use Only. Not for use in diagnostic procedures.