

Anti-DENND2D Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A12572-2-Fluoro488

About DENND2D

Enables guanyl-nucleotide exchange factor activity. Located in cytosol and nucleoplasm.

Overview

Product Name	Anti-DENND2D Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9H6A0

Technical Details

Immunogen	E.coli-derived human DENND2D recombinant protein (Position: M1-K471).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DENND2D Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.