

## Anti-DENND2D Antibody Picoband® HRP Conjugated

Catalog Number: A12572-2-HRP

### About DENND2D

Enables guanyl-nucleotide exchange factor activity. Located in cytosol and nucleoplasm.

### Overview

|                      |   |
|----------------------|---|
| Product Name         | Anti-DENND2D Antibody Picoband® HRP Conjugated                                      |
| Reactive Species     | Human   |
| Clonality            | Polyclonal  |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.    |
| Host                 | Rabbit  |
| Uniprot ID           | Q9H6A0  |

### Technical Details

|                     |  |
|---------------------|--|
| Immunogen           | E.coli-derived human DENND2D recombinant protein (Position: M1-K471).                              |
| Form                | Liquid   |
| Concentration       | 0.5 mg/mL  |
| Purification        | Immunogen affinity purified.   |
| Conjugate           | HRP  |
| Suggested Dilutions | The intended application should be selected according to the customer's experimental requirements. |

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.