

## Anti-MRP-S12 Antibody

Catalog Number: A12643S12


### Overview

Product Name	Anti-MRP-S12 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-MRP-S12 Antibody catalog # A12643S12. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O15235

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 69-114 of Human MRP-S12.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

## Anti-MRP-S12 Antibody (A12643S12) Images



Western blot (WB) analysis of MRP-S12 (K59) pAb at 1:1000 dilution Lane1:AML-12 whole cell lysate(40ug) Lane2:H9C2 whole cell lysate(40ug) Lane3:HCT116 whole cell lysate(40ug) Lane4:L02 whole cell lysate(40ug) Lane5:HeLa whole cell lysate(40ug)

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MRP-S12 Antibody

For Research Use Only. Not for use in diagnostic procedures.