

## Anti-TNNI3K FPGT Antibody

Catalog Number: A12778


### Overview

Product Name	Anti-TNNI3K FPGT Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TNNI3K FPGT Antibody catalog # A12778. Tested in ELISA, IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O14772


### Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human TNNI3K. AA range:301-350
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 IHC 1:100-1:300 ELISA 1:10000 IF 1:50-200

## Anti-TNNI3K FPGT Antibody (A12778) Images



Western blot analysis of Mouse-kidney lysis using TNNI3K antibody. Antibody was diluted at 1:500. Secondary antibody was diluted at 1:20000



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using TNNI3K Antibody. The picture on the right is blocked with the synthesized peptide.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TNNI3K FPGT Antibody

For Research Use Only. Not for use in diagnostic procedures.