

Anti-SKOR2 Antibody

Catalog Number: A12823

About SKOR2

SKOR2 (SKI family transcriptional corepressor 2), also named FUSSEL18, is a transcriptional repressor that is highly expressed in the central nervous system (PMID: 16200078). SKOR2 was detected at 30 and 60 kD in the cerebellum. It is identified as a novel transcriptional regulator in Purkinje cells that acts upstream of Shh during cerebellum development (PMID: 21937600). Skor2 was detected in Purkinje cells in the human fetal cerebellum (PMID: 34140698).

Overview

Product Name	Anti-SKOR2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-SKOR2 Antibody catalog # A12823. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q2VWA4

Technical Details

Immunogen	Peptide
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunofluorescence, 1:50-400 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SKOR2 Antibody

For Research Use Only. Not for use in diagnostic procedures.