

Anti-Zygote arrest protein 1 ZAR1 Antibody

Catalog Number: A12936-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-Zygote arrest protein 1 ZAR1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Zygote arrest protein 1 ZAR1 Antibody catalog # A12936-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q86SH2

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 220-270 of Human ZAR1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-Zygote arrest protein 1 ZAR1 Antibody (A12936-1) Images



Western blot (WB) analysis of ZAR1 polyclonal antibody at 1:500 dilution Lane1:The Testis tissue lysate of Mouse(40ug) Lane2:The Testis tissue lysate of Rat(40ug) Lane3:K562 whole cell lysate(40ug) Lane4:Hela whole cell lysate(40ug) Lane5:SGC7901 whole cell lysate(40ug) Lane6:The Kidney tissue lysate of Mouse(40ug) Lane7:The Kidney tissue lysate of Rat(40ug)

Anti-Zygote arrest protein 1 ZAR1 Antibody

For Research Use Only. Not for use in diagnostic procedures.