

Anti-Epiphycan EPYC Antibody

Catalog Number: A12937

About EPYC

Essential for female fertility. May play a role in the oocyte-to-embryo transition By similarity.

Wu X., Nat. Genet. 33:187-191(2003).

Wu X., Biol. Reprod. 69:861-867(2003).

Overview

Product Name	Anti-Epiphycan EPYC Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Epiphycan EPYC Antibody catalog # A12937. Tested in ELISA, WB applications. This antibody reacts with Human, Mouse.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q99645

Technical Details

Immunogen	Synthesized peptide derived from human EPYC protein.
Predicted Reactive Species	Chimpanzee, Drosophila, Macaque
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	EPYC antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the

optimal dilution ratio for your samples.
Some PubMed article(s) citing the expression level of this target are as follows:
Boster Bio's internal QC testing used:
WB 1:500-2000 ELISA 1:5000-20000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Epiphycan EPYC Antibody