

Anti-RDH14 Antibody

Catalog Number: A12970

About RDH14

Enables all-trans-retinol dehydrogenase (NADP+) activity. Involved in osteoblast differentiation. Located in cytosol; endoplasmic reticulum; and nucleoplasm.

Overview

| | |
|----------------------|--|
| Product Name | Anti-RDH14 Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-RDH14 Antibody catalog # A12970. Tested in WB, IHC applications. This antibody reacts with Human, Mouse, Rat. |
| Application | IHC, WB |
| Clonality | Polyclonal |
| Formulation | 500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | 12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing. |
| Host | Rabbit |
| Uniprot ID | Q9HBH5 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human RDH14 recombinant protein (Position: V3-V332). |
| Form | Liquid |
| Concentration | 500 ug/ml |
| Purification | Immunogen affinity purified. |
| Suggested Dilutions | Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RDH14 Antibody

For Research Use Only. Not for use in diagnostic procedures.