

Anti-MSTO1 Antibody Picoband® Cy3 Conjugated

Catalog Number: A13028-1-Cy3

About MSTO1

Involved in mitochondrion distribution; mitochondrion organization; and positive regulation of mitochondrial fusion. Located in mitochondrial outer membrane. Is active in cytosol.

Overview

Product Name	Anti-MSTO1 Antibody Picoband® Cy3 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9BUK6

Technical Details

Immunogen	E.coli-derived human MSTO1 recombinant protein (Position: E44-Q518). Human MSTO1 shares 78.1% amino acid (aa) sequence identity with mouse MSTO1.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-MSTO1 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.