

Anti-MSTO1 Antibody Picoband® PE Conjugated

Catalog Number: A13028-1-PE

About MSTO1

Involved in mitochondrion distribution; mitochondrion organization; and positive regulation of mitochondrial fusion. Located in mitochondrial outer membrane. Is active in cytosol.

Overview

| | |
|----------------------|--|
| Product Name | Anti-MSTO1 Antibody Picoband® PE Conjugated |
| Reactive Species | Human |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q9BUK6 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human MSTO1 recombinant protein (Position: E44-Q518). Human MSTO1 shares 78.1% amino acid (aa) sequence identity with mouse MSTO1. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-MSTO1 Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.