

Anti-ELAC1 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A13042-2-Fluoro550

About ELAC1

Enables tRNA-specific ribonuclease activity. Involved in rescue of stalled ribosome and tRNA 3'-end processing. Located in cytosol and nucleoplasm.

Overview

Product Name	Anti-ELAC1 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9H777

Technical Details

Immunogen	E.coli-derived human ELAC1 recombinant protein (Position: F35-D353).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ELAC1 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.