

Anti-ELAC1 Antibody Picoband® PE Conjugated

Catalog Number: A13042-2-PE

About ELAC1

Enables tRNA-specific ribonuclease activity. Involved in rescue of stalled ribosome and tRNA 3'-end processing. Located in cytosol and nucleoplasm.

Overview

Product Name	Anti-ELAC1 Antibody Picoband® PE Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9H777

Technical Details

Immunogen	E.coli-derived human ELAC1 recombinant protein (Position: F35-D353).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-ELAC1 Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.