

Anti-SGT2/SGTB Antibody Picoband® Biotin Conjugated

Catalog Number: A13056-1-Biotin

About SGTB

Predicted to enable molecular adaptor activity. Predicted to be involved in positive regulation of chaperone-mediated protein folding; posttranslational protein targeting to endoplasmic reticulum membrane; and ubiquitin-dependent ERAD pathway. Predicted to be part of TRC complex. Predicted to be active in membrane.

Overview

Product Name	Anti-SGT2/SGTB Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q96EQ0

Technical Details

Immunogen	E.coli-derived human SGT2/SGTB recombinant protein (Position: Q20-A268).
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SGT2/SGTB Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.