

Anti-Glypican 2/GPC2 Antibody Picoband® Cy3 Conjugated

Catalog Number: A13075-2-Cy3

About GPC2

Predicted to be involved in several processes, including positive regulation of neuron projection development; regulation of protein localization to membrane; and smoothed signaling pathway. Located in endoplasmic reticulum.

Overview

| | |
|----------------------|--|
| Product Name | Anti-Glypican 2/GPC2 Antibody Picoband® Cy3 Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q8N158 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E.coli-derived human Glypican 2/GPC2 recombinant protein (Position: E25-D498). Human Glypican 2/GPC2 shares 83.3% amino acid (aa) sequence identity with both mouse and rat Glypican 2/GPC2. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-Glypican 2/GPC2 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.