

## Anti-Glypican 2/GPC2 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A13075-2-Fluoro647

### About GPC2

Predicted to be involved in several processes, including positive regulation of neuron projection development; regulation of protein localization to membrane; and smoothed signaling pathway. Located in endoplasmic reticulum.

### Overview

Product Name	Anti-Glypican 2/GPC2 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8N158

### Technical Details

Immunogen	E.coli-derived human Glypican 2/GPC2 recombinant protein (Position: E25-D498). Human Glypican 2/GPC2 shares 83.3% amino acid (aa) sequence identity with both mouse and rat Glypican 2/GPC2.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

**Anti-Glypican 2/GPC2 Antibody - Fluoro647**

For Research Use Only. Not for use in diagnostic procedures.