

## Anti-Glypican 2/GPC2 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A13075-2-Fluoro647

### About GPC2

Predicted to be involved in several processes, including positive regulation of neuron projection development; regulation of protein localization to membrane; and smoothed signaling pathway. Located in endoplasmic reticulum.

### Overview

Product Name	Anti-Glypican 2/GPC2 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8N158

### Technical Details

Immunogen	E.coli-derived human Glypican 2/GPC2 recombinant protein (Position: E25-D498). Human Glypican 2/GPC2 shares 83.3% amino acid (aa) sequence identity with both mouse and rat Glypican 2/GPC2.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

**Anti-Glypican 2/GPC2 Antibody - Fluoro647**

For Research Use Only. Not for use in diagnostic procedures.