

Anti-Glypican 2/GPC2 Antibody Picoband® HRP Conjugated

Catalog Number: A13075-2-HRP

About GPC2

Predicted to be involved in several processes, including positive regulation of neuron projection development; regulation of protein localization to membrane; and smoothed signaling pathway. Located in endoplasmic reticulum.

Overview

Product Name	Anti-Glypican 2/GPC2 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8N158

Technical Details

Immunogen	E.coli-derived human Glypican 2/GPC2 recombinant protein (Position: E25-D498). Human Glypican 2/GPC2 shares 83.3% amino acid (aa) sequence identity with both mouse and rat Glypican 2/GPC2.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Glypican 2/GPC2 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.