

Anti-TTLL10 Antibody Picoband® Cy3 Conjugated

Catalog Number: A13080-1-Cy3

About TTLL10

Predicted to enable protein-glycine ligase activity, elongating. Predicted to be involved in protein polyglycylation. Predicted to be located in axoneme and microtubule cytoskeleton.

Overview

Product Name	Anti-TTLL10 Antibody Picoband® Cy3 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6ZVT0

Technical Details

Immunogen	E.coli-derived human TTLL10 recombinant protein (Position: S205-P658). Human TTLL10 shares 65.4% and 64.9% amino acid (aa) sequence identity with mouse and rat TTLL10, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TTLL10 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.