

## Anti-TTLL10 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A13080-1-Fluoro594

### About TTLL10

Predicted to enable protein-glycine ligase activity, elongating. Predicted to be involved in protein polyglycylation. Predicted to be located in axoneme and microtubule cytoskeleton.

### Overview

Product Name	Anti-TTLL10 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6ZVT0

### Technical Details

Immunogen	E.coli-derived human TTLL10 recombinant protein (Position: S205-P658). Human TTLL10 shares 65.4% and 64.9% amino acid (aa) sequence identity with mouse and rat TTLL10, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

**Anti-TTLL10 Antibody - Fluoro594**

**For Research Use Only. Not for use in diagnostic procedures.**