

Anti-beta-1,3-Gal-T4 B3GALT4 Antibody

Catalog Number: A13136T4


Overview

Product Name	Anti-beta-1,3-Gal-T4 B3GALT4 Antibody		
Reactive Species	Human, Mouse, Rat		
Description	Boster Bio Anti-beta-1,3-Gal-T4 B3GALT4 Antibody catalog # A13136T4. Tested in ELISA, WB applications. This antibody reacts with Human.		
Application	Application		ELISA, WB
Clonality	Polyclonal		
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation. For conjugation methods that are not compatible with the stabilizing component of the formulation, a carrier-free antibody format is required.		
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to 6 months. Avoid repeated freeze-thaw cycles.		
Host	Rabbit		
Uniprot ID	O96024		

Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human B3GALT4. AA range:181-230
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Anti-beta-1,3-Gal-T4 B3GALT4 Antibody (A13136T4) Images



Western blotting validation for Anti-beta-1

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-beta-1,3-Gal-T4 B3GALT4 Antibody

For Research Use Only. Not for use in diagnostic procedures.