

Anti-Thioredoxin-like protein 4B TXNL4B Antibody

Catalog Number: A13142


Overview

Product Name	Anti-Thioredoxin-like protein 4B TXNL4B Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Thioredoxin-like protein 4B TXNL4B Antibody catalog # A13142. Tested in WB,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NX01


Technical Details

Immunogen	Recombinant full length Human TXNL4B.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IF: 1:50-1:200

Anti-Thioredoxin-like protein 4B TXNL4B Antibody (A13142) Images



Western blot analysis of extracts of mouse liver, using TXNL4B Antibody.



Immunofluorescence analysis of A549 cell using TXNL4B antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Thioredoxin-like protein 4B TXNL4B Antibody

For Research Use Only. Not for use in diagnostic procedures.