

Anti-RUFY2 Antibody Picoband® APC Conjugated

Catalog Number: A13178-1-APC

About RUFY2

RUN and FYVE domain containing 2 (RUFY2) is a protein that in humans is encoded by the RUFY2 gene. Enables SH3 domain binding activity. Predicted to be involved in regulation of endocytosis. Located in cytoplasm.

Overview

| | |
|----------------------|--|
| Product Name | Anti-RUFY2 Antibody Picoband® APC Conjugated |
| Reactive Species | Human |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q8WXA3 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human RUFY2 recombinant protein (Position: E239-K551). |
| Cross Reactivity | No cross-reactivity with other proteins |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RUFY2 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.