

Anti-RUFY2 Antibody Picoband® HRP Conjugated

Catalog Number: A13178-1-HRP

About RUFY2

RUN and FYVE domain containing 2 (RUFY2) is a protein that in humans is encoded by the RUFY2 gene. Enables SH3 domain binding activity. Predicted to be involved in regulation of endocytosis. Located in cytoplasm.

Overview

Product Name	Anti-RUFY2 Antibody Picoband® HRP Conjugated
Reactive Species	Human
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8WXA3

Technical Details

Immunogen	E.coli-derived human RUFY2 recombinant protein (Position: E239-K551).
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RUFY2 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.