

Anti-TPPP3 Antibody Picoband® APC Conjugated

Catalog Number: A13192-1-APC

About TPPP3

Enables tubulin binding activity. Involved in decidualization and microtubule bundle formation. Located in microtubule bundle and perinuclear region of cytoplasm.

Overview

Product Name	Anti-TPPP3 Antibody Picoband® APC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9BW30

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human TPPP3.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TPPP3 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.