

## Anti-TPPP3 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A13192-1-Fluoro594

### About TPPP3

Enables tubulin binding activity. Involved in decidualization and microtubule bundle formation. Located in microtubule bundle and perinuclear region of cytoplasm.

### Overview

Product Name	Anti-TPPP3 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9BW30

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human TPPP3.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TPPP3 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.