

Anti-SYT13 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A13202-2-Fluoro647

About SYT13

This gene encodes a member of the large synaptotagmin protein family. Family members have an extracellular N-terminal transmembrane domain and a cytoplasmic C terminus with two tandem C2 domains (C2A and C2B). Synaptotogmin family members can form homo- and heteromeric complexes with each other. They also have different biochemical properties and developmental profiles, and patterns of tissue distribution. Synaptotagmins function as membrane traffickers in multicellular organisms. Two alternatively spliced transcript variants that encode different protein isoforms have been described for this gene.

Overview

| | |
|----------------------|--|
| Product Name | Anti-SYT13 Antibody Picoband® Fluoro647 Conjugated |
| Reactive Species | Mouse, Rat |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q7L8C5 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human SYT13 recombinant protein (Position: S19-L426). |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SYT13 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.