

Anti-MEIOB Antibody Picoband® APC Conjugated

Catalog Number: A13211-2-APC

About MEIOB

Predicted to enable chromatin binding activity; single-stranded DNA 3'-5' exodeoxyribonuclease activity; and single-stranded DNA binding activity. Predicted to be involved in double-strand break repair via homologous recombination; fertilization; and meiotic nuclear division. Predicted to be located in cytoplasm. Implicated in spermatogenic failure 22.

Overview

Product Name	Anti-MEIOB Antibody Picoband® APC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8N635

Technical Details

Immunogen	E.coli-derived human MEIOB recombinant protein (Position: I140-V442). Human MEIOB shares 77.9% and 76.1% amino acid (aa) sequence identity with mouse and rat MEIOB, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MEIOB Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.