

Anti-MEIOB Antibody Picoband® Biotin Conjugated

Catalog Number: A13211-2-Biotin

About MEIOB

Predicted to enable chromatin binding activity; single-stranded DNA 3'-5' exodeoxyribonuclease activity; and single-stranded DNA binding activity. Predicted to be involved in double-strand break repair via homologous recombination; fertilization; and meiotic nuclear division. Predicted to be located in cytoplasm. Implicated in spermatogenic failure 22.

Overview

Product Name	Anti-MEIOB Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8N635

Technical Details

Immunogen	E.coli-derived human MEIOB recombinant protein (Position: I140-V442). Human MEIOB shares 77.9% and 76.1% amino acid (aa) sequence identity with mouse and rat MEIOB, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MEIOB Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.