

## Anti-C18orf22/RBFA Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A13257-1-Fluoro550

### About RBFA

Predicted to be involved in rRNA processing. Predicted to be located in mitochondrion.

### Overview

Product Name	Anti-C18orf22/RBFA Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8N0V3

### Technical Details

Immunogen	E.coli-derived human C18orf22/RBFA recombinant protein (Position: L15-D319).
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-C18orf22/RBFA Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.