

## Anti-Phospholipase ABHD3 ABHD3 Antibody

Catalog Number: A13293

### About ABHD3

Has no peptidase activity.

Quesada V., m. Biophys. Res. Commun. 314:54-62(2004).

Puente X.S., Genome Res. 14:609-622(2004).

### Overview

Product Name	Anti-Phospholipase ABHD3 ABHD3 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Phospholipase ABHD3 ABHD3 Antibody catalog # A13293. Tested in ELISA, WB applications. This antibody reacts with Human, Mouse.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8WU67

### Technical Details

Immunogen	Synthesized peptide derived from human ABHD3 protein.
Predicted Reactive Species	Chimpanzee, Drosophila, Macaque
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ABHD3 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the

optimal dilution ratio for your samples.  
Some PubMed article(s) citing the expression level of this target are as follows:  
Boster Bio's internal QC testing used:  
WB 1:500-2000 ELISA 1:5000-20000

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospholipase ABHD3 ABHD3 Antibody